


# Introduction to Cyclo Olefin Polymer (COP)

## - Key Properties Update -

S. Suzuki, T. Katayama, K. Arai, T. Sawaguchi, ZEON CORPORATION

### 1. What is Cyclo Olefin Polymer ?

◆ COP polymers (ZEONEX®, ZEONOR®) were commercialized in 1990 and have found increasing use in pharmaceutical syringes and vials due to COP's unique benefits vs glass and other plastics.



Product Grade	Water Absorption (%)	Light transmittance (%)	Glass Transition Temp (°C)	Elongation at break (%)
ZEONEX® 5000	<0.01	92	69	120
ZEONOR® 1020R	<0.01	92	102	90
ZEONEX® 690R	<0.01	92	136	20
ZEONEX® 790R	<0.01	92	163	10

Passes US/EU/Japan Pharmacopoeia, ISO 10993, DMF listed

### 2. Benefits of COP for pre-filled syringe

Syringes made of ZEONEX® offer :


- ✓ High transparency  
- Easy to inspect the drug
- ✓ High break resistance  
- Improved drop tolerance
- ✓ Low impurities  
- Very low residual metals
- ✓ Low E / L  
- Very low elution
- ✓ Sterilization  
- Gamma, EB, Steam
- ✓ High moisture barrier  
- Long-term drug storage
- ✓ No delamination  
- No flakes and particles
- ✓ Low adsorption  
- Keep drug activity
- ✓ Drug compatibility  
- Acids, alkalis and alcohols
- ✓ Silicone oil free  
- No aggregation

Typical devices where ZEONEX® is used:

- Pre-Filled syringes
- Pre-Filled cartridges
- Vials and bottles for long-term storage of biologics
- IV and Total Parenteral Nutrition (TPN) bags
- Bio-reactors
- High pressure injection syringes (Needle-free, viscous drug, and large dosage injection)

### 3. Tough Material

In case of syringe shape, ZEONEX® has high strength compared with glass.



Test speed of 1.0 mm/min


### 4. Safe Material

- High Purity
- Residual metals less than 0.02 ppm
- Contains no lubricants/process aids
- Low outgas

Residual metals of ZEONEX® 690R						
Al	V	Zr	Mg	Ti	Pd	
<0.02	<0.02	<0.02	<0.02	<0.02	<0.02	
Cr	W	Fe	Ni	Zn	Cd	
<0.02	<0.02	<0.02	<0.02	<0.02	<0.02	

Unit: ppm  
Measures by ICP-MS analysis

Very low residual metals (vs. other plastics)




Measurement method: ICP-MS


Analysis of metal species : Al, Cu, Fe, Mg, Pb, Ti, Zn, Zr, Li

### 5. Low Extractables / Leachables

◆ Non-polar components



◆ Polar components



Non-polar and polar components are not detected from ZEONEX® syringes.

Test conditions : Storage temp. : 40°C, Storage time : 3 months (approximately equal to storage of 3 years at 5°C)

### 6. Sterilize with Gamma, EB, Steam

ZEONEX® is minimally influenced by exposure to standard sterilization methods.


Properties	Unit	Initial	Sterilization #2			
			Steam	EOG	Gamma 25 kGy	EB 20 kGy
Light Transmittance *1	%	91	91	91	87	88
Yellow Index (ΔYI) *1	-	0.3	0.8	0.3	4.1	5.8
HAZE #1	-	0.1	0.6	0.2	0.1	0.2
Tensile Strength	MPa	68	74	68	67	67

\*1 Test piece: injection-molded plate (3 mm in thickness)

\*2 Sterilization conditions

Steam: 121 deg.C, 20 min. EOG: 50° C, 6 hours, EOG conc. 600 mg/l


Color shift occurs after irradiation sterilization but quickly recovers



### 7. Low Moisture Pass - Thru


- Excellent moisture barrier → prolong drug shelf life
- 3x better gas barrier vs polypropylene (PP)

Moisture barrier



Test conditions  
Measurement method : JIS K7129A

O<sub>2</sub>, N<sub>2</sub> barrier




Test conditions  
Measurement method : JIS K7126A, at 23 °C


### 8. No Delamination - No flakes and particles -

- Delamination most commonly occurs at the tip of glass syringe where process heat history exposure is the highest.
- Particles are observed with glass syringes.
- No particles are observed in ZEONEX® syringe.

Syringe made of ZEONEX®



Glass syringe




Test conditions

3.0 % Citric Acid, pH : 10.0, Storage temp. : 80 °C, Storage time : 28 days


More details

Identification of particles



Measurement method : EDX (Energy Dispersive X-ray)

Glass surface profile after testing



The particles may be generated from the cavities observed on the glass surface.

- Identification of the particles
- Capture of the particles.